



Arise



Redefining patient handling

Supportive, adaptable and extremely flexible

Arise

The Arise was specifically designed for quick, simple and safe transfers. Unlike using conventional passive lifting devices, the Arise gives support by using a range of active slings. This encourages the patient to participate during the lifting process.



This active style of transfer is proven to keep patients partially weight bearing for longer. Slings are simple to fit, provide excellent access and support for the following key tasks:

- Toileting
- Standing
- Walking/Rehab
- Quick transfers

The Arise is a highly adaptable lifting device and is capable of transferring a patient from a low seated position, to fully extended with its unrivalled range of movement.



The Arise is packed with useful features such as:

Adjustable knee pad

The knee pad is height adjustable by 6"/152mm for added flexibility. This ensures the patient feels comfortable and fully supported during the lifting process. The pad is moulded for additional support and there is a leg restraint for those patients needing additional reassurance.



Removable foot tray

Removing the foot tray transforms the Arise into a walking aid which can facilitate rehabilitation and training. With the foot tray removed you can fully raise a patient to a standing position. You can then remove the Arise and allow the patient to walk.



The Arise was subjected to a rigorous usability study before its release. This type of user feedback has ensured that the Arise is robust, usable and flexible enough to suit the needs of the patient, carer and the environment in which it operates.

Dedicated support handle



The handle is carefully positioned to ensure that optimum comfort and safety can be achieved whilst the patient is being transferred.

5

The Arise was designed with the 5 key principles of moving and handling in mind. This direct correlation ensures Oxford's lift designs are simple, safe and above all comfortable to use.

This view is shared by our clinical consultant teams, who have assisted us in conducting detailed usability surveys to ensure that the Arise meets the true needs of the patient, the carer and the environment in which it operates.

- ① Get close to the load
- ② Use a wide stable base
- ③ Ensure a comfortable firm grasp
- ④ Keep spine close to neutral
- ⑤ Make sure movement is smooth

Specification	Imperial	Metric
Safe Working Load	375 lbs	170 kgs
Maximum Overall Length	53.9 inches	1370 mm
Minimum Overall Length	51.1 inches	1300 mm
Maximum Overall Height	67.7 inches	1720 mm
Minimum Overall Height	44.0 inches	1120 mm
Maximum Height to Attachment Point	66.9 inches	1700 mm
Minimum Height to Attachment Point	40.1 inches	1020 mm
Turning Radius	55.9 inches	1420 mm
Legs Open - External Width	40.5 inches	1030 mm
Legs Open - Internal Width	35.0 inches	890 mm
Legs Closed - External Width	24.8 inches	630 mm
Legs Closed - Internal Width	20.0 inches	510 mm
Overall Height of Legs	4.7 inches	120 mm
Ground Clearance	1.4 inches	35 mm
Front Twin Castors	4.0 inches	100 mm
Rear Braked Castors	4.0 inches	100 mm

Weights	Imperial	Metric
Mast, Base & Boom Assembly	108.9 lbs	49.5 kgs
Power Pack	6.6 lbs	3.0 kgs
Total	115.5 lbs	52.5 kgs
Base Assembly	42.9 lbs	19.5 kgs
Mast & Boom Assembly	51.7 lbs	23.5 kgs

Specifications may vary due to manufacturing tolerances

Electrical Specifications

Battery 1 X 12 volt rechargeable sealed lead acid type
 Battery capacity 3.2 Ampere hours
 Charger rated input 230Vac 50/60Hz
 Charger rated output 27.4/29.0 VDC@0.8A

Warranty

The Arise comes with a 2 year warranty as standard.
 Please refer to your user manual for further information.

Electrical Shock Protection

Charger - class II
 Lift - internal power source

Degree of shock protection
 Charger - Type B
 Lift - Type B

Intended operating environment : >+5° <+40°
 Outside this environment functionality and safety may be compromised



Joerns Healthcare Limited
 High Street, Wollaston, Stourbridge, West Midlands DY8 4PS England
 (T) +44 (0)1384 446622 (F) +44 (0)1384 446601
 www.joerns.co.uk • email: info@joerns.co.uk